

## Appendix C

# Interferogram and Profile analysis

### C.1 Python script: Profile calculation for CSI

The extraction of a surface profile is done following the method of CSI. The whole process is discussed in Sec. 3.2. The **basic operations** are the following:

```

1  """
2  Load interferogram files and create a stack
3  """
4  """
5  ...
6  """
7  """
8  Load calibration file
9  """
10 """
11 ...
12 """
13
14 def gaus(x,a,x0,sigma,c):
15     return a*np.exp(-(x-x0)**2/(2*sigma**2))+c
16
17 """
18 PER PIXEL ANALYSIS
19 """
20 for i in range(rows):
21     for j in range(cols):
22         # extract interferogram in vertical axis (normal to the surface)
23         z = image_stack[:,i,j]
24         zm = z - np.median(z)
25         za = np.abs(zm)
26         # take the "envelope"; sigma is selected as a good trade-off
27         # with visual criteria.
28         zg = gaussian_filter1d(za, sigma=30)
29
30         # make a first guess for the gaussian fitting of the "envelope"
31         nmax = np.argmax(zg)
32         amp = zg[nmax]
33
34         try:
35             # statistically 60 pixels is a good guess for the specific
36             # white light source used
37             sigma = 60
38             c = np.median(zg)
39             guess = [amp,nmax,sigma,c]
40             popt, _ = curve_fit(gaus,n,zg,p0=guess)
41             peak = popt[1]

```

```

40         # the surface elevation is the interpolated value of the
        calibration array with input
41         # the calculated peak of the best-fit .
42         profile[i,j] = curve(peak)
43     except:
44         profile[i,j] = None
45
46     """
47 Save profile
48     """
49     """
50     ...
51     """

```

## C.2 Python script: Profile calculation for SRWLI

This process is presented in Sec. E. The **basic operations** are the following:

```

1
2     """
3 Load images of interference , reference , sample and background.
4 Load array containing the calibrated values for the decomposed
  wavelengths.
5     """
6
7 # compute wavenumber array
8 k = 2*np.pi/lamda
9
10 # smooth the following signals before using them to isolate the
    interference cosine term
11 Is_g = gaussian_filter1d(Is , sigma=3, axis=1)
12 Ir_g = gaussian_filter1d(Ir , sigma=3, axis=1)
13 Ib_g = gaussian_filter1d(Ib , sigma=3, axis=1)
14
15 # calculate interference cosine term --> 'Ipr'
16 eps = np.finfo(np.float64).eps
17 nom = 2.0*np.sqrt(Is_g)*np.sqrt(Ir_g)
18 nom[np.where(nom==0)] = eps
19 Ipr = (Itot-Is_g-Ir_g-Ib_g)/nom
20
21 # EXTRA PRE-PROCESSING
22 # extract waviness --> 'Iprmodg2'
23 Ipr_g = gaussian_filter1d(Ipr , sigma=11, axis=1)
24 Ipr_mod = Ipr - Ipr_g
25 Ipr_modg = gaussian_filter1d(Ipr_mod , sigma=3, axis=1)
26 Iprmodg2 = Ipr_modg - np.median(Ipr_modg)
27
28 # Linear resampling (factor=1) of the wavenumber samples
29 num_resample = cols
30 if num_resample % 2 == 0:
31     num_resample -= 1
32 k_new = np.linspace(k[0] , k[-1], num_resample)
33
34 # calculate frequency range and samples
35 dk = abs(k_new[0] - k_new[1])
36 Fn = 1/(2*dk)
37 Nfreqs_f = num_resample
38 f_axis = np.linspace(0, Fn, Nfreqs_f//2+1)
39
40 # initialize profile array
41 z = np.zeros(rows)

```

```

42 """
43 Iterate procedure for each row of the image
44 """
45 for i in range(rows):
46     # resampled waviness signal
47     f_rsmpl = interp1d(k, Iprmodg2[i], kind='cubic')
48     Ipr_new = f_rsmpl(k_new)
49
50     # FT calculation
51     fourier = np.fft.fft(Ipr_new, n=Nfreqs_f)
52     fourier = np.abs(fourier)
53     fourier = np.fft.fftshift(fourier)
54     fourier = fourier[Nfreqs_f//2:]
55
56     # Upsampling of fourier signal to "increase" resolution
57     Nf_times = 20
58     f_axisnew = np.linspace(f_axis[0], f_axis[-1], Nfreqs_f*Nf_times)
59     f_rsmpl2 = interp1d(f_axis, fourier, kind='cubic')
60     fourier_new = f_rsmpl2(f_axisnew)
61
62     # smooth fourier signal to eradicate noise for peak detection
63     sigma_f = Nf_times*1 #Nf_times/4
64     fourier_g = gaussian_filter(fourier_new, sigma=sigma_f)
65     argmax_idx = np.argmax(fourier_g)
66     # calculate surface elevation
67     z[i] = np.pi*f_axisnew[argmax_idx]/um
68

```

### C.3 Python script: Piezo actuator calibration

This process is explained thoroughly in Sec. 3.6. The **basic operations** are the following:

```

1 """
2 Load interferogram files and create a stack
3 """
4 ...
5 """
6
7
8 """
9 PHASE CALCULATION FROM EACH INTERFEROGRAM
10 """
11 phases[0] = calc_phase(image)
12
13 for j in range(1, listlen):
14     image = ((io.imread(file_list[j])).astype(float))
15     # odd dimensions of image —> Frequency 0 of Fourier axes is exactly
16     # on the center of the image
17     image = im2odd_dim(image)
18     phases[j] = calc_phase(image)
19
20 """
21 CALCULATE DISPLACEMENT FROM PHASE SHIFTS
22 """
23 uwphases = np.unwrap(phases)
24 dp = np.diff(uwphases)
25
26 # Phase difference converted to displacement of fringes.
27 # In michelson interferometers, displacement of the sample is equal to
28 # half the optical path difference

```

```

28 lamda = float(input('Enter the peak wavelength of the source (nm)...\n'))
    ) # in nm
29 displacements = dp/(2*np.pi)*lamda/2 # in nm
30 displacements = np.concatenate([[0],displacements]) # zero as the first
    position
31
32 incremental = np.cumsum(displacements)
33 grad = np.gradient(incremental)
34 grad = reject_outliers(grad, m = 2.)
35 dp_mean = np.mean(grad)
36
37 incremental = - incremental if dp_mean < 0 else incremental
38
39 """
40 Save Calibration file
41 """
42 ...
43 """

```

## C.4 Python script: Profile post-processing

This process is explained in Sec. 3.3. It requires user visual feedback and for this reason a small GUI has been built. User events signal callback functions that do the basic operations of the GUI. The **callback functions** are the following:

```

1 """
2 Create Matplotlib window with widgets for buttons, figures, axes, span
    selectors and ROI selectors.
3 """
4 """
5 ...
6 """
7
8 def line_select_callback(eclick, erelease):
9     global ROI
10
11     if toggle_selector_RS.active:
12         x1, y1 = int(round(eclick.xdata)), int(round(eclick.ydata))
13         x2, y2 = int(round(erelease.xdata)), int(round(erelease.ydata))
14         ROI = [(x1,y1), (x2,y2)]
15
16 def enable_ROIselect(event):
17     global ts_id
18
19     del ROI[:]
20
21     toggle_selector_RS.set_active(True)
22     toggle_selector_RS.set_visible(True)
23     ts_id = fig.canvas.mpl_connect('key_press_event', toggle_selector)
24
25 def toggle_selector(event):
26     if event.key=='enter' and toggle_selector_RS.active:
27         # print(' RectangleSelector deactivated. ')
28         textbox.set_val('ROI selected')
29         toggle_selector_RS.set_visible(False)
30         toggle_selector_RS.update()
31         toggle_selector_RS.set_active(False)
32         fig.canvas.mpl_disconnect(ts_id)
33
34 def onselect(xmin, xmax):
35     global flag_save, thres, unrotated

```

```

36     thres[thres<xmin] = xmin
37     thres[thres>xmax] = xmax
38     ax1.cla()
39     ax1.hist(thres.ravel(), bins=255)
40     ax1.set_title('Press left mouse button and drag to select the wanted
41                  histogram range')
42     ax2.imshow(thres, cmap=plt.cm.viridis)
43     ax2.set_title('Thresholded profile')
44     unrotated = 'thres'
45     textbox.set_val('Low-high thresholds set as [%d,%d]'%(xmin,xmax))
46     cbar2.set_clim(vmin=np.amin(thres),vmax=np.amax(thres))
47     cbar2.draw_all()
48     flag_save = 1
49
50 def crop(event):
51     global thres, flag_save
52
53     (x1,y1), (x2,y2) = ROI
54     thres = thres[y1:y2,x1:x2]
55
56     ax1.cla()
57     ax1.hist(thres.ravel(), bins=255)
58     ax1.set_title('Press left mouse button and drag to select the wanted
59                  histogram range')
60     ax2.imshow(thres)
61     ax2.set_title('Cropped profile')
62     cbar2.set_clim(vmin=np.amin(thres),vmax=np.amax(thres))
63     cbar2.draw_all()
64     textbox.set_val('Image cropped')
65     flag_save = 3
66
67 def rotate_image(event):
68     global flat_profile, flag_save, unrotated
69
70     (x1,y1), (x2,y2) = ROI
71     if unrotated == 'thres':
72         flat_profile = rotate(thres, [y1, y2], [x1, x2])
73         im = ax1.imshow(thres)
74         ax1.axis('off')
75         ax1.set_title('Original/Thresholded')
76         cbar1 = fig.colorbar(im, ax=ax1)
77         cbar1.set_clim(vmin=np.amin(thres),vmax=np.amax(thres))
78         cbar1.draw_all()
79     elif unrotated == 'raw':
80         flat_profile = rotate(prof, [y1, y2], [x1, x2])
81         im = ax1.imshow(prof)
82         ax1.axis('off')
83         ax1.set_title('Original/Raw')
84         cbar1 = fig.colorbar(im, ax=ax1)
85         cbar1.set_clim(vmin=np.amin(prof),vmax=np.amax(prof))
86         cbar1.draw_all()
87     ax2.imshow(flat_profile)
88     ax2.set_title('Rotated/Horizontalized')
89     cbar2.set_clim(vmin=np.amin(flat_profile),vmax=np.amax(flat_profile))
90     cbar2.draw_all()
91     textbox.set_val('Image rotated')
92     flag_save = 2
93
94 def compute_parameters(event):
95     global fig2

```

```

96     (x1,y1), (x2,y2) = ROI
97     area = flat_profile[y1:y2,x1:x2]
98     # area = thres[y1:y2,x1:x2]
99     # area = prof[y1:y2,x1:x2]
100
101     # Mean plane
102     M = np.mean(area)
103
104     # Distances of the suface from the mean plane
105     abs_diff = np.abs(area - M)
106
107     # Flatness deviation
108     # flt =
109
110     # Maximum peak height of the surface
111     Sp = np.amax(area-M)
112
113     # Maximum valley depth of the surface
114     Sv = np.amax(M-area)
115
116     # Maximum height of the surface
117     Sz = Sp + Sv
118
119     # Average roughness ( Mean difference around mean plane )
120     Sa = np.mean(abs_diff)
121
122     # Roughness (Standard deviation of difference(smoothed surface ,
123     surface))
124     filt = gaussian_filter(area, 11)
125     diff = filt - area
126     Sr = np.std(diff)
127
128     # Root mean square roughness
129     Sq = np.sqrt(np.mean(abs_diff**2))
130
131     # Skewness (Degree of symmetry of the surface heights around the
132     mean plane)
133     Ssk = 1/Sq**3*(np.mean(abs_diff**3))
134
135     # Kurtosis (Presence of lack of inordinately high peaks or deep
136     valleys)
137     Sku = 1/Sq**4*np.mean(abs_diff**4)
138
139     col_labels = ['Surface\nParameter', 'Value', 'Unit', 'Description']
140     table_data = np.array([[ 'M', '%.0f'%M, 'nm', 'Mean plane'], [ 'Sp', '%.0f'
141     '%Sp, 'nm', 'Maximum peak height of the surface'],
142     [ 'Sv', '%.0f'%Sv, 'nm', 'Maximum valley depth of
143     the surface'], [ 'Sz', '%.0f'%Sz, 'nm', 'Maximum height of the surface'
144     ],
145     [ 'Sa', '%.0f'%Sa, 'nm', 'Roughness: Mean
146     difference around mean plane'],
147     [ 'Sr', '%.0f'%Sr, 'nm', 'Roughness: Standard
148     deviation around surface waviness'],
149     [ 'Sq', '%.0f'%Sq, 'nm', 'Roughness: Root mean
150     square difference around mean plane'],
151     [ 'Ssk', '%.4f'%Ssk, '-', 'Skewness: Degree of
152     symmetry of the surface heights around the mean plane.\n'+
153     'Ssk>0: peaks predominance, Ssk
154     <0: valleys predominance'],
155     [ 'Sku', '%.4f'%Sku, '-', 'Kurtosis: Presence (
156     Sku>3) of lack (Sku<3) of '+
157     'inordinately high peaks or deep
158     valleys']]

```

```

146
147 # MATPLOTLIB TABLE (NOT SELECTABLE TEXT)
148
149 fig2 , ax3 =plt.subplots(1,1)
150 figManager = plt.get_current_fig_manager()
151 figManager.window.showMaximized()
152 fig2.canvas.set_window_title('Areal Height Parameters')
153 ax3.axis('tight')
154 ax3.axis('off')
155 tparam = ax3.table(cellText=table_data , colLabels=col_labels ,
156                   colWidths=[0.1,0.1,0.1,0.75] ,
157                           cellLoc='center' , colLoc='center' , loc='center')
158 tparam.auto_set_font_size(False)
159 tparam.set_fontsize(12)
160 tparam.scale(1,1.6)
161
162 def average_pplot(event):
163     (x1,y1) , (x2,y2) = ROI
164     area = flat_profile[y1:y2,x1:x2]
165
166     ax1.cla()
167     if direction == 'horiz profile':
168         avg = np.average(area , axis=0)
169         ax1.plot(avg)
170         ax1.set_title('Average horizontal profile of ROI')
171     else:
172         avg = np.average(area , axis=1)
173         ax1.plot(avg)
174         ax1.set_title('Average vertical profile of ROI')

```

## C.5 Python script: Radius of Curvature assessment (Curved profile)

This process is explained in Sec. 3.5. The **basic operations** are the following:

```

1 """
2 Load profile image
3 """
4 """
5 ...
6 """
7
8 # values in um
9 # calibrated pixel size range
10 px_range = [9.8 , 10.2]
11 px_mean = 10.0
12 px_uc = 0.2
13 sigma_z = 0.028
14
15 sel = '0'
16
17 while sel != '1' and sel != '2':
18     print('NOMINAL AND MEASURED RADIUS OF CURVATURE COMPARISON\n\n')
19     print('Enter \'1\' if nominal ROC is known...\n')
20     print('Enter \'2\' if refractive index and focal distance of lens
21         are known...\n')
22
23     sel = input()
24
25 if sel == '1':
26     ROC = float(input('Enter ROC (mm)\n'))

```

```

26     f = float(input('Enter nominal focal distance (mm)\n'))
27 elif sel == '2':
28     n_lens = float(input('Enter nominal refractive index of the lens...\n'))
29     f = float(input('Enter nominal focal distance (mm)\n'))
30     ROC = abs(f*(n_lens-1.0))
31
32     """
33     BEST-FIT CALCULATION OF 1) 2D PROFILE WITH A PARABOLOID
34                     2) 1D PROFILE PASSING THROUGH CENTER WITH A PARABOLA
35     FOR UPPER AND LOWER LIMIT OF PIXEL SIZE
36     """
37
38     for i in range(2):
39         px = px_range[i]
40
41         fringearea_gf = gaussian_filter(fringearea, sigma=1)
42
43         if f < 0:
44             [r0], [c0] = np.where(fringearea_gf == np.amin(fringearea_gf))
45             Z = fringearea - fringearea[r0, c0]
46             radius_c = ROC*1000
47         elif f > 0:
48             [r0], [c0] = np.where(fringearea_gf == np.amax(fringearea_gf))
49             Z = fringearea - fringearea[r0, c0]
50             radius_c = - ROC*1000
51
52
53         X = np.arange(c)*px
54         Y = np.arange(r)*px
55         X, Y = np.meshgrid(X, Y)
56
57         guess = [radius_c, px*r0, px*c0, 0.0]
58
59         Z = Z/1000 # now in um
60
61         po, pc = curve_fit(paraboloid, (X,Y), Z.ravel(), p0=guess)
62         paraboloid_fit = paraboloid((X,Y), *po).reshape(np.shape(Z))
63         perr = np.sqrt(np.diag(pc))
64
65         po_list.append(po)
66         pc_list.append(pc)
67         perr_list.append(perr)
68
69         if f < 0:
70             [r0], [c0] = np.where(fringearea_gf == np.amin(fringearea_gf))
71             z = fringearea[r0]-np.amin(fringearea[r0])
72             radius_c = ROC*1000
73         elif f > 0:
74             [r0], [c0] = np.where(fringearea_gf == np.amax(fringearea_gf))
75             z = fringearea[r0]-np.amax(fringearea[r0])
76             radius_c = - ROC*1000
77
78         guess2 = [radius_c, px*c0]
79
80         z = z/1000 # now in um
81
82         x = np.arange(c)*px
83         s_z = sigma_z*np.ones(c)
84         po2, pc2 = curve_fit(parabola, x, z, p0=guess2, sigma=s_z,
85                             absolute_sigma=False)
86         perr2 = np.sqrt(np.diag(pc2))

```



```

87     po2_list.append(po2)
88     pc2_list.append(pc2)
89     perr2_list.append(perr2)
90
91     """
92     STATISTICS
93     """
94
95     residuals = Z.ravel() - paraboloid((X,Y), *po)
96     ss_res = np.sum(residuals**2)
97     ss_tot = np.sum((Z.ravel()-np.mean(Z.ravel()))**2)
98     r_squared = 1 - (ss_res / ss_tot)
99
100     r_squared_list.append(r_squared)
101
102     """
103     RESULTS: ROC AND UNCERTAINTY IN MEASUREMENT
104     """
105     # values in mm
106     error_pxmin = perr_list[0][0]/1000
107     error_pxmax = perr_list[1][0]/1000
108     ROCmeas_pxmin = abs(po_list[0][0])/1000
109     ROCmeas_pxmax = abs(po_list[1][0])/1000
110     minROC = ROCmeas_pxmin - error_pxmin
111     maxROC = ROCmeas_pxmax + error_pxmax
112     ROC_meas = (maxROC + minROC)/2
113     ROC_error = (maxROC - minROC)/2

```



## Appendix D

# Integrated Graphical User Interface (GUI)

For the sake of operational simplicity of the whole procedure from the interferogram capturing to the profile analysis, a minimalistic GUI was developed (Fig:D.1, D.2). The requirements of the preparation for measurement are such that many procedures has to be done in parallel. For this reason, the code of this GUI was built strongly based on the "multiprocessing" package of Python. To run the GUI, simply run the "WLtopo\_mp.py" script from the Windows command prompt (some IDEs may be problematic with running multi-processing modules in full functionality).




FIGURE D.1: GUI main window: operations and graph planes for results.




FIGURE D.2: GUI camera configuration window: all basic options needed for the camera to run.

**—FUNCTIONS—****CAMERA CONFIGURATION TAB**

Before taking a measurement, many things must be done in order, so as to achieve the best possible measurement. For starters, the user connects the camera, initializes it and configures it according to the "config.ini" file. A copy of it is stored in an editable text box in the "camera configuration tab". The user can change every single camera option of each section. There are three sections: "*Default*", which uses the maximum possible image size, "*Image Sequence Configuration*", which uses 16-bit pixel depth and internal software trigger for capturing each image and "*Live View Configuration*", which does not use any trigger but has 8-bit pixel depth for displaying purposes. Any other camera option is irrelevant to the occasion used and freely changes by just overwriting and saving the corresponding values.

**LIVE VIEW**

On the left of the "Operations" tab, the "Live View" window displays continuously the selected region of interest of the sample. Press "Live!" and "Stop Live feed" to start and stop respectively the live feed.

**PREPERATION FRAME***Select ROI*

All operations here have the purpose of setting up the sample holder of the interferometer in the best possible position in the x,y,z axes, with the most important being the axis parallel to the optical one.

- Click on the "Select ROI" button.
- Click on the "Live View" WIndow on the wanted top left corner of the ROI and release on the bottom right one. A box selection has been made.
- Modify the selection by clicking and releasing around the box.
- When ready, press enter.

*Display intensity across horizontal line*

In some cases when seeking interference fringes, it might be useful to have an auxiliary live plot of just a line across the image to search for deviations in the constant light intensity that may indicate the existence of fringes. In that case:

- Press "Display intensity across horizontal line". Light intensity vs horizontal image axis plot is displayed on the right window/plane in synchronization with the live view feed on the left window.
- Press the same button to stop the plotting.

*Oscillate Piezo*

WLSI method requires the existence of a well shaped interference pattern per sample point (signal per pixel) from which the argmax will be easily determined, inferring thus the height of this sample point. So to avoid marginal situations where the argmax is located very close to the beginning or the end of the "per pixel signal", an auxiliary function can be used.

The user can select multiple Points of Interest (POI):

- Click on the "Select POI" button.
- Click on the "Live View" window on the wanted POI. Press enter.
- If more points of interest are needed, select each one and press enter.
- When every POI is selected, to finish the process click on the "OK" button next to the "Select POI" button.
- Press "Live!" to start the live feed.
- Press "Oscillate Piezo" to start the back and forth oscillation of the piezo. In each movement of the piezo,

## MEASUREMENT FRAME

*Piezo scan – Capture interferograms*

Following preparation, for the WLSI measurement and calibration of the piezo movements, a piezo scanning must be done for each case:

- Type a measurement/calibration folder name (which will contain the interferograms) on the corresponding text field.
- Select "Measurement" or "Calibration" from the radio buttons.
- Click on "Piezo scan – Capture interferograms" button to execute the image acquisition. [Note: The calibration process requires the proper use of the optical arrangement (Coherent light source, optical flat/mirror in reference arm)]

*Calibrate*

Before the first measurement, a calibration must be done regarding the piezo movements, that is the steps and real distance made correspondence. [Note: The calibration process requires the proper use of the optical arrangement (Coherent light source, optical flat/mirror in sample arm)]. After the acquisition of the images as described above, the user will:

- Click on the "Calibrate" button. When the process is done, a text file containing the piezo displacements between each scan will be saved on the root folder with name "Mapping\_Steps\_Displacement\_<piezo steps>\_<file increment number>.txt"

*Analyze*

The final step to the whole process for obtaining a 2D profile of a specimen is the analysis of the measurement interferograms:

- Click on the "Analyze" button.
- Select the interferograms image set from the emerging prompt window.
- When the profile is calculated, it will be saved on the root folder with name "raw\_profile\_<interferograms image set folder name>" in text and tiff format.
- A new window will emerge for *post processing* of the profile. These functions are *image cropping*, *histogram thresholding*, *average line profiling* and *rotation of the profile plane parallel to the ground*.

**–PARALLELIZATION–**

The workload specifically for live continuous image display and synchronized feature plotting is beyond the capabilities of usual single-threaded script/modules. In this way, each new operation starts only when the previous in order has finished. To overcome this, a multi-processing/multi-threading approach was implemented which allows the parallel execution of at least two different operations.

Two extra processes are created when running the GUI script "WLtopo\_mp.py". More specifically, there exist a process for *GUI design and display* which actually creates a child process for *camera manipulation* purposes and another for *piezo oscillation* purposes. In order to have synchronized camera displaying, piezo movement and feature plotting, these three processes must communicate with each other. Every time some data are needed to be plotted either on "Live View" window or in the adjacent plane (feature plotting), the camera manipulation process transfer these data to a shared queue. Then, the GUI process reads the data from the queue and displays/plots them accordingly. The whole operation happens in a datum per datum manner to keep a synchronized data "production" and "consumption".